

Submitting dental framework designs

Introduction

This document gives guidance on submitting framework designs from open CAD systems to Renishaw's central manufacturing. It assumes that the design has been created and saved as a '.stl' file to a central location on your pc and that the pc has internet access as well as Renishaw inciseCAD installed and ready for use.

Note: It is only possible to manufacture crowns and bridges in cobalt chrome from '.stl' files.

Instructions

Open Renishaw inciseCAD 📜 1.



Click Windows Start > All Programs > Renishaw PLC > inciseCAD > inciseCAD


Click on New Prescription 2.



Click on Add Framework 3.



4. Define the framework.



Select the framework type - 'Crown / Bridge'.

Select the material - 'Cobalt Chrome'.


Note that the 'Thickness' and 'Cement Space' options have no effect on the design.

- 5. Click OK.
- Select Import CAD Data 6.





7. Browse to wherever the '.stl' design file is saved and click **Open**. You will see the framework design in the main view.



In the prescription explorer tree, you will see that prescription is given the 'green light' indicating 8. that it is ready to be submitted ${\color{orange} {\mathbb{Q}}}$

Ensure the prescription is highlighted and click on **Submit To Central Manufacturing**



Click 'Yes' to the dialog that opens and the design is submitted.

Prescription status

You can get live updates on the prescription status by using the 'Prescription Management' function in inciseCAD.

Click on the 'Prescription Management' tab and in the explorer pane, highlight 'Recent Central Submissions'. A list of all submitted prescriptions is displayed. Click on the prescription in the prescription tree and a dialog will appear showing the status of the prescription.



If the prescription has been dispatched you can track it by clicking on Courier Link.

