

Metrology leader launches new Polish website

World leading metrology company Renishaw, has launched a new Polish language website to support the growth of its customer base in Poland.

The site contains details on the full range of Renishaw products, including CMM probe systems, touch probes and laser tool setters for CNC machine tools, linear and rotary encoders for position feedback, calibration systems for machine performance analysis, styli for touch probes, Raman spectroscopy systems for spectral analysis and 3D digitising products for reverse engineering.

Visitors to the site will also be able to download a range of Polish language brochures, data sheets, news articles and technical articles.

Renishaw Sp. z o.o. was formed in Warsaw in 2002, and National Sales Manager, Tomek Rzyško has welcomed the new web site. "There are few international suppliers in the metalworking industry with Polish websites, and this new website highlights Renishaw's commitment to educate and support the Polish market."

There are now Renishaw web sites in nineteen different languages. As well as the major European languages, sites are also available in Korean, Traditional Mandarin, Simplified Mandarin, Czech, Dutch, Russian, Swedish and Slovenian.

The Renishaw Group has 33 offices operating in 27 countries, with a total of 1,865 employees.